

Main Features

	Agil Pro	Agil
Speech Guard	•	•
Binaural Noise Management	•	-
Binaural Processing (compression)	•	-
Binaural Synchronization (automatics)	•	-
Binaural Coordination (push-buttons)	•	•
Binaural Dynamic Feedback Cancellation 2 (DFC2)	•	•
Multi-Band Adaptive Directionality	•	•
TriState Noise Management	•	•
Power Bass (streaming)	•	•
Music Widening (streaming)	•	•
Life Learning	•	•



Fitting Options (Genie 2012.1)

	Agil Pro	Agil
nEARcom Cordless enabled	•	•
In-Situ audiometry	•	•
VAC Rationale	•	•
NAL-NL1, NAL-NL2, DSL v 5.0a m[i/o]	•	•
Identities	5	5
Fitting Bands	10	10
Activity Analyzer	•	•

Additional Features

	Agil Pro	Agil
Fitting Bandwidth*	10 kHz	10 kHz
3 RITE Receivers (S, M and P)	•	•
Front Focus	•	•

Connectivity

	Agil Pro	Agil
ConnectLine Enabled	•	•
DAI and FM	•	•



*Bandwidth accessible for gain adjustments during fitting

miniRITE



RITE



CUSTOM



HS/LP



ITC/HS Power (Omni or Dir)*



ITC



CIC/MIC Power*



CIC/MIC*

* Not wireless

miniBTE



Hook



Corda² thin tube

BTE13



Hook



Corda² thin tube

BTE POWER



Hook

SPEAKER AND COUPLING OPTIONS



Open Dome



Plus Dome



Power Dome



LiteTip



Micro Mold



Power Micro Mold



Power Mold

CORDA² OPTIONS




Open Dome




Plus Dome




Custom Mold



STANDARD



MEDIUM



POWER

■ Custom Mold - Power Dome

□ Open Dome - Plus Dome

■ HS/LP

□ CIC/MIC Power

■ Earmold

□ Corda²

RITE & BTE INSTRUMENTS

90 Chroma Beige

94 Terracotta

93 Chestnut Brown

44 Silver

91 Silver Grey

92 Steel Grey

63 Diamond Black

74 Samoa Blue

▶ miniRITE style only

CUSTOM INSTRUMENTS

01 Pink

02 Tan

03 Medium Brown

04 Dark Brown