



Features	Oticon Agil Pro	Oticon Agil
Software	Genie 2010.1	Genie 2010.1
ITE Styles: Power Cic through Low Profile	Yes	Yes
BTE Styles: MiniRITE/RITE/BTE 13	Yes	Yes
3 RITE Speakers (S,M,P)	Yes	Yes
Binaural Connectivity	Binaural Processing	Binaural program and volume change
Spatial Sound 2 with Binaural Compression Binaural Noise Management	Yes	No
Speech Guard	Yes	Yes
Streamer/ConnectLine	Yes	Yes
Power Bass	Yes	Yes
Music Widening	Yes	Yes
Binaural DFC2 Directionality	Yes: Automatic Adaptive Multiband (4) 3 modes	Yes: Automatic Adaptive Multiband (4) 3 modes
Noise Management	TriState Noise Management with Voicefinder	TriState Noise Management with Voicefinder
Front Focus	RITE/BTE styles	RITE/BTE Styles
My Voice	Yes	No
Bandwidth	10,000Hz	10,000Hz
Adjustable Channels	10	10
Identities	5	5
Fitting Rationales	Voice Aligned Compression NAL-NL1;DSL v5.0	Voice Aligned Compression NAL-NL1;DSL v5.0
Life Learning (VC learning)	Yes	Yes
Analyzer (Data Logging)	Yes	Yes
Number of programs	4	4
AutoPhone/T-Coil except miniRITE/cic and mic	Yes	Yes
FM Compatible	RITE/BTE	RITE/BTE
Conda Compatible	BTE 13	BTE 13
nEARcom	All models except cic	All models except cic