



906 58 611 00/1-09

Oticon ♦ Dual



Great sound
never looked better



People First

We believe that it takes more than technology and audiology to create the best hearing instruments. That's why we put the individual needs and wishes of people with hearing loss first in our development of new hearing care solutions.

- 1) Bandwidth 10 kHz in Dual M9 model, 8 kHz bandwidth in Dual M7 and M5 models.
- 2) Multiband adaptive directionality available in Dual M9 and M7 models. Single band directionality available in M5 model.
- 3) Tristate Noise Management in Dual M9 and M7 models feature. Two state noise management in M5 model.

www.oticonusa.com



Inside Dual Mini



A beauty on the outside A sensation on the inside

Feel more like yourself again

Extended Bandwidth

Dual captures a wider range of sound than other hearing devices, particularly in the high frequencies, giving you more speech information cues e.g. from words with consonants such as “s”, “f” and “th”. This extended bandwidth³⁾ gives you a unique advantage as it improves ease of listening and speech understanding.

Relax and join in the conversation

Adaptive directionality

Directionality in a hearing device helps you to focus on speech in noisy situations. Dual supports the world’s most advanced directional microphone system²⁾, which means that you can hear people talking even when other noise sources get in the way. Dual can be set to automatically reduce noise coming from behind you as well as from the sides. Even if you move or turn your head, Dual will continue to reduce noise so that you can concentrate fully on what is being said. Enjoy superior speech understanding even in challenging situations. With Dual, you can relax and enjoy yourself at restaurants and dinner parties again.

Experience speech without distractions

Superior noise management

Dual’s noise management system³⁾ works in perfect harmony with its directional system. This ensures that any sudden, unpleasant or dominant noises will be automatically reduced. With Dual, you will always be free to focus on what is being said without having to make a great effort in noisy surroundings.

Wind noise protection

In addition, Dual features a special wind filter that protects the directional microphones, giving you a more comfortable experience even in windy conditions.

Eliminate Feedback

Dynamic Feedback Cancellation² (DFC²)

Dual contains a new, advanced and extremely fast reacting feedback cancellation system that detects and eliminates potential feedback before it’s even audible.

Adjust easily

Automatic Adaptation Manager

All hearing instruments require a little getting used to. Dual’s Automatic Adaptation Manager is specifically designed to help you to adapt more easily over time. By gradually increasing the level of amplification, automatically until it reaches the target set by your hearing care professional, Dual will help you to adjust comfortably and automatically.

Benefit from greater reliability

Nanocoating

Dual also features nanocoating. This special anti-moisture treatment is designed to protect your hearing devices and improve reliability.

Look forward to phone conversations

AutoPhone

AutoPhone will make phone calls much easier. When you put your phone to your ear, Dual will switch automatically to a customized program that is tuned for your listening needs. Whether you use a home or office phone, or a cell phone, you will quickly access that call without any extra handling. After the call, Dual will automatically switch back to your normal setting.

Efficient battery performance

Optimized battery life

Dual was designed to give you the optimal combination of advanced performance, discreet size and battery life. Dual gives you approximately 105 hours per battery so you will have efficient battery performance.

Enjoy a hearing device just for you

Choice of colors and performance options

Dual was designed specifically for those who want the most advanced technology available in the smallest package. Available in your choice of 14 colors, Dual Mini instruments are the most discreet possible but are fully equipped with the sound processing power available in today’s more advanced hearing instruments. Experience who you can be again. It’s not just what Dual can do for you. But what you can do with Dual.

