


EPOQ IN DETAIL

ONE COMPLETE PRODUCT FAMILY, UNLIMITED FITTING OPPORTUNITIES



EPOQ XW

THE ULTIMATE HEARING SOLUTION

Featuring spatial sound for the most realistic hearing possible

Binaural Processing XW 10 KHz

Two Epoqs work in unison constantly processing and exchanging data to secure optimal settings for the levels and compression resulting in best possible interaural time and level differences — a pivotal element of spatial sound perception.

EPOQ W

HIGH END HEARING SOLUTION

Featuring binaural synchronization and wireless connectivity

Binaural Synchronization W 10 KHz

The high-speed wireless technology intelligently synchronizes automatics, noise reduction and directionality to preserve speech intelligibility in all listening situations.

EPOQ V

HIGH PERFORMANCE HEARING CARE

Featuring high performance hearing care with wireless connectivity at a more affordable price

Binaural Coordination V 8 KHz

What represents the limit of wireless technologies in other products today forms the basis of Epoq. Full coordination of volume and program shifts ensure a baseline for a balanced sound picture. The simplest approach to wireless.

	EPOQ XW	EPOQ W	EPOQ V
SPATIAL SOUND			
Extended Bandwidth	10 KHz	10 KHz	8 KHz
Binaural Processing	●	○	○
Binaural Synchronization	●	●	○
Binaural Coordination	●	●	●
Binaural DFC 2	●	○	○
DFC 2	●	●	●
True Dynamics	●	●	○
Front Focus	●	●	●
My Voice	●	○	○
RATIONALE			
VAC	●	●	○
NAL/DSL	○	○	●
LIFE LEARNING			
Compression & volume control	●	●	Only VC
CONNECTIVITY			
Streamer	●	●	●

● Yes ○ No



NEW COLOR PALETTE

Standard Line



Chroma Beige



Chestnut Brown



Silver Grey



Steel Grey

Black Line



Chroma Beige



Espresso



White Silver



Steel Grey

White Line



Gold Dust



Ice Blue



Orchid



Silver Grey

Cool Colors



Blue



Red



Purple



Silver



Baby Pink



Baby Blue

RITE



Micro Mold



Plus Dome



Open Dome

Power RITE



Power Dome



Power Mold

THE ULTIMATE ACCESSORY

Connect conveniently and directly to a range of communications and entertainment devices



Remote Control

The Streamer is now a fully functional Remote, allowing for convenient control of volume and program shifts.

Wireless Phone Headset

Connect the Streamer to a cell phone or landline phone using a standard Bluetooth protocol. In effect, you have turned a pair of Epoqs into a hands-free, wireless headset.

MP3 Headphones

Input any MP3 player directly into the Streamer to turn a pair of Epoqs into a wireless headset.

Home Entertainment Link

Connect the Streamer to a variety of entertainment devices either via Bluetooth or directly:

- TV
- Radio
- Stereo System
- PC or laptop

NEW Silicone Skins

Customize your Streamer with new protective colored silicone skins

- Dark Grey
- White
- Blue
- Orange
- Green

