

Features & Overview

| Technical Features: | Hit Pro | Hit |
|---|-------------|-------|
| Bandwidth | 8 kHz | 8 kHz |
| RISE Platform* | ✓ | ✓ |
| Automatic Directionality (Surround & Split) | ✓ | ✓ |
| Adaptive Directionality | Single band | |
| Noise Management (modulation) | ✓ | ✓ |
| Dynamic Feedback Cancellation 2 | ✓ | ✓ |
| Open Ear Acoustics | ✓ | ✓ |
| Automatic Adaptation Manager | ✓ | |
| AutoPhone program | ✓ | ✓ |
| Front Focus | ✓ | ✓ |
| Wind noise protection | ✓ | ✓ |
| Fitting (Genie 2009.1): | | |
| NAL-NL1 & DSL v5.0a m[i/o] | ✓ | ✓ |
| Fitting Bands | 6 | 4 |
| No. of user programs | 4 | 4 |
| nEARcom cordless enabled (RITE/BTE only) | ✓ | ✓ |
| Identities | 1 | 1 |
| Activity Analyzer | ✓ | ✓ |
| Live View | ✓ | ✓ |
| More: | | |
| T-coil | ✓ | ✓ |
| Corda ² thin tube solution | ✓ | ✓ |
| DAI & FM | ✓ | ✓ |
| Battery Low Warning | ✓ | ✓ |
| Sound indicators for program shifts | ✓ | ✓ |
| On-set delay and jingle | ✓ | ✓ |
| Mute/Standby mode | ✓ | ✓ |

* Featuring superior sound quality and automatics using outcome-based decision making and parallel processing.

Oticon ♦ Hit Pro /Hit

Styles

RITE 312



RITE

BTE 312



Earhook



Corda² thin tube

BTE 13



Earhook



Corda² thin tube

BTE 13 Power



Earhook

Custom Instruments



CIC/MIC



ITC



HS/LP

Colors

Chroma Beige



Steel Grey



Silver Grey



Chestnut Brown



Purple



Blue



Baby Blue



Red



Baby Pink



Silver

