

# Oticon ♦ Epoq



Small **SIZE**  
superb **SOUND**  
smart **CONNECTIVITY**

**oticon**  
PEOPLE FIRST




# Get the HOTTEST technology the MAX sound quality and an instant CONNECTION - with EPOQ!

With some of the hottest technology on planet Earth, Epoq connects you to your friends - and your gadgets - with the ultimate sound quality imaginable.

With Epoq you get detail, depth and balance. So you can pinpoint the sounds you want to hear, communicate with all your friends and get closer than ever to that ultimate goal: a natural listening experience.

With smart, discreet devices on your ears and the wireless sound Streamer around your neck, Epoq becomes your 21st-century gateway to a brand new world of sound.



**TechWizardry**

*Where design, sound quality and connectivity are concerned, Epoq takes hearing aids to a new level. Never before have such powerful sound processing technologies been packed into such a small and smart device.*



## Outdoors or indoors be FREE to COMMUNICATE and stay TOTALLY in the PICTURE

Some days, chatting on your cell phone can be pretty tough, but with Epoq you can be completely confident. You don't have to hold your phone to your hearing aid or worry about interference, because the Bluetooth-enabled Streamer keeps you wirelessly in the picture.


---

### **TechWizardry**

*The Epoq Streamer is a Bluetooth device that acts as a gateway to external sound sources. It turns the Epoq hearing instruments into wireless ear-phones with super sound quality!*







# Super-clear SPEECH mega-sweet MUSIC and the power to EXPRESS yourself

When it comes to sound quality, Epoq is invincible. As the most hi-tech hearing optimizer ever to be produced, it makes speech sound super clear and music mega sweet . Whether you're studying at school or jamming with friends, Epoq gives you the power to express yourself.

**TechWizardry**

*With its extreme bandwidth, Epoq captures difficult-to-hear high-frequency sounds. These softer sounds add depth and detail to the sound picture. Speech and music become clearer, and everyday sounds become easier to localize.*






# When you're on the MOVE and in the MOOD to GROOVE turn your Epoqs into FUNKY EARPHONES

With Epoq you can get lost in the music, any time and anywhere!  
Connect the Streamer to your MP3 and it'll turn your tiny  
Epoq instruments into the funkiest wireless earphones. With  
awesome sound you'll get more out of rock, pop or rap tracks,  
whether you're at home, on a train - or up a tree!

**TechWizardry**

*The Streamer makes it very easy to manage electronic devices. If the phone should ring while you're listening to music, it can be answered at the mere touch of a button. The Streamer even has its own Volume Control.*

A person is shown from the chest up, wearing a black headset with a microphone. They are looking towards the right side of the frame, presumably at a screen. The background is blurred, showing what appears to be a gaming setup with red and black components. The overall lighting is soft and focused on the person.

# CHILL out join the GAME and kiss clumsy cables GOODBYE

Whether you're chilling out to music or gaming with friends, you want great sound quality all the time. With the Streamer you can forget about wires, too, because it pipes all the digital details straight into both ears.


---

## **TechWizardry**

*As a companion device, the Streamer sends the signal wirelessly, directly to both ears. This helps to strengthen the signal to noise ratio and keeps interference to a minimum.*





A person is shown from the chest up, wearing a white neck-worn device (Streamer) and a small case (Epoq). The person is wearing a red shirt and green pants. The background is a plain white wall.

## Look COOL feel CONFIDENT and impress all your FRIENDS

While the small Epoq casing is a neat and discreet shape, the neck-worn Streamer is a real showstopper. This hi-tech duo will make you look cool and feel confident, all day long. When your friends find out how the Streamer allows hands-free gadget management, they'll be totally envious!

**TechWizardry**

*Epoq's compact design, ultra-thin sound tube and tiny receiver in-the-ear make it a highly cosmetic solution. And the Streamer is as elegant as any MP3 player.*





## Be TRENDY be CHOOSY be YOURSELF!



You can get Epoqs that sit snugly inside your ear, and Epoqs that sit firmly behind it. These come in a range of cool metallic colors designed to suit your personality. These devices are so easy to wear that you're likely to forget all about them.



**TechWizardry**

*Sophisticated manufacturing processes create Epoq's cool colors and sturdy finishes, while top quality components ensure long-term reliability.*



**People First** We believe that it takes more than technology and audiology to create the best hearing instruments. That's why we put the individual needs and wishes of people with hearing loss first in our development of new hearing care solutions.