

		Epoq Power			Vigo Connect Power		Hit Power	
		XW	W	V	Vigo Pro Connect	Vigo Connect	Hit Pro	Hit
BTE Power	MPO: 127 dB Pk	●	●	●	●	●	●	●
	Gain: 61 dB Pk							
AUDIOLOGICAL FEATURES								
	Bandwidth	10 kHz	10 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz
	Binaural Compression System	●	○	○	○	○	○	○
	Binaural Synchronization	●	●	○	○	○	○	○
	Binaural Coordination	●	●	●	●	●	○	○
	Binaural DFC2	●	○	○	○	○	○	○
	DFC2	●	●	●	●	●	●	●
	True Dynamics	●	●	○	○	○	○	○
	Directionality	Multiband Adaptive	Multiband Adaptive	Multiband Adaptive	Multiband Adaptive	Adaptive	Adaptive	Automatic, Fixed
	Noise management	TriState/ Voice Finder	TriState/ Voice Finder	TriState/ Voice Finder	TriState/ Voice Finder	Modulation based	Modulation based	Modulation based
	Identities	5	5	3	3	3	1	1
	Fitting bands	10	10	8	8	6	6	4
	T-coil	●	●	●	●	●	●	●
	Automatic Adaptation Manager	●	●	●	●	●	●	○
	Battery Size/ Life	13/215 hrs.	13/215 hrs.	13/215 hrs.	13/215 hrs.	13/215 hrs.	13/215 hrs.	13/215 hrs.
RATIONALES								
	Voice Aligned Compression (VAC)	●	●	○	○	○	○	○
	NAL NL-1/DSL v5.0a	○	○	●	●	●	●	●
LIFE LEARNING								
	Compression & Volume Control	●	●	Only VC	Only VC	○	○	○
WIRELESS CONNECTIVITY								
	Connectivity	●	●	●	●	●	○	○
OPTIONS								
	DAI	●	●	●	●	●	●	●
	FM shoe	R12	R12	R12	R12	R12	R12	R12
	FM adapter	FM9	FM9	FM9	FM9	FM9	FM9	FM9

YES = ● ● ●

NO = ○ ○ ○

RE 2 CC Coupler