

EPOQ RITE POWER

VIGO CONNECT RITE POWER

RE2 CC COUPLER	XW	W	V	VIGO PRO CONNECT	VIGO CONNECT
MPO: 124dB PK	●	●	●	●	●
Gain: 55dB PK	●	●	●	●	●
Battery Size/Life	312/100 hrs	312/100 hrs	312/100 hrs	312/100 hrs	312/100 hrs

AUDIOLOGICAL FEATURES

Extended Bandwidth	10 KHz	10 KHz	8 KHz	8 KHz	8 KHz
Spatial Sound Compression System	●	○	○	○	○
Binaural Synchronization	●	●	○	○	○
Binaural Coordination	●	●	●	●	●
Binaural DFC 2	●	○	○	○	○
DFC 2	●	●	●	●	●
True Dynamics*	●	●	○	○	○
Directionality	MultiBand Adaptive	MultiBand Adaptive	MultiBand Adaptive	MultiBand Adaptive	Adaptive
Noise Management	TriState/VoiceFinder	TriState/VoiceFinder	TriState/VoiceFinder	TriState/VoiceFinder	Modulation
DAI	●	●	●	●	●

RATIONALE

Voice Aligned Compression (VAC)**	●	●	○	○	○
NAL/DSL	○	○	●	●	●

LIFE LEARNING

Compression & Volume Control	●	●	Only VC	Only VC	○
------------------------------	---	---	---------	---------	---

WIRELESS CONNECTIVITY†

Streamer	●	●	●	●	●
----------	---	---	---	---	---

● Yes ○ No

● Yes ○ No

* True Dynamics combines a slow acting and fast acting compression system. This protects listening quality at all sound levels, maximizing audibility and speech understanding while minimizing fatigue.

** VAC: Oticon's proprietary fitting rationale—curvilinear compression strategy with multiple knee-points. Compression reaches the highest point in the soft/moderate input area and becomes increasingly linear as input levels increase. Improves speech understanding especially at high input levels.

† FM can be connected via the Streamer or the T-coil with a neckloop.