



Oticon Rise Hearing Instruments

	Epoq XW	Epoq W	Epoq V	Vigo Pro Connect	Vigo Connect	Hit Pro	Hit
ITE Styles: CIC through Low Profile	Yes, Except CIC /MIC	Yes	Yes	Yes	Yes	Yes	Yes
BTE Styles: RITE/Power RITE /BTE/Power BTE	Yes	Yes	Yes	Yes	Yes	Yes, except RITE power	Yes, except RITE power
Binaural Connectivity	Binaural Processing	Binaural Synchronization	Binaural Coordination	Binaural Coordination	Binaural Coordination	No	No
Streamer/ConnectLine	Yes	Yes	Yes	Yes	Yes	No	No
Directionality	Automatic MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Automatic MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Automatic MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Automatic MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Automatic Single Band Adaptive Directionality with 3 modes: surround, split, full	Automatic Single Band Adaptive Directionality with 2 modes: surround, split	Automatic Single Band Fixed Directionality with 2 modes: surround, split
Noise Management	TriState Noise Management with Voicefinder	TriState Noise Management with Voicefinder	TriState Noise Management with Voicefinder	TriState Noise Management with Voicefinder	Speech-weighted modulation based	Speech-weighted modulation based	Speech-weighted modulation based
Bandwidth (Hz)	10,000	10,000	8,000	8,000	8,000	8,000	8,000
Adjustable Channels	10	10	8	8	6	6	4
DFC2	Yes (binaural)	Yes	Yes	Yes	Yes	Yes	Yes
Fitting Rationales	Voice Aligned Compression (VAC)	Voice Aligned Compression (VAC)	NAL-NL1; DSL v5.0	NAL-NL1; DSL v5.0	NAL-NL1; DSL v5.0	NAL-NL1; DSL v5.0	NAL-NL1; DSL v5.0
Identities	5	5	3	3	3	1	1
Front Focus	BTE & RITE styles	BTE & RITE styles	BTE & RITE styles	BTE & RITE styles	BTE & RITE styles	BTE & RITE styles	BTE & RITE styles
My Voice	Yes	No	No	No	No	No	No
Life Learning (VC Learning)	Life Learning	Life Learning	VC Learning	VC Learning	No	No	No
Analyzer (Data Logging)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Number of Programs	4	4 (1 for CIC)	4 (1 for CIC)	4 (1 for CIC)	4 (1 for CIC)	4 (1 for CIC)	4 (1 for CIC)
AutoPhone/T-coil (except CIC/MIC)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
FM Compatibility	BTE Only	BTE Only	BTE Only	BTE Only	BTE Only	BTE Only	BTE Only
Corda Compatible (Open Fit/Thin Tube)	BTE 13 hook (except power)	BTE 13 hook (except power)	BTE 13 hook (except power)	BTE 312 / 13 hook (except power)	BTE 312 / 13 hook (except power)	BTE 312 / 13 hook (except power)	BTE 312 / 13 hook (except power)
nEARcom	All models except CIC	All models except CIC	All models except CIC	All models except CIC	All models except CIC	RITE/BTE only	RITE/BTE only



Oticon ♦ Dual Oticon Dual Family

	Dual Connect XW	Dual Connect W	Dual Connect V	Dual Mini 9	Dual Mini 7	Dual Mini 5
Spatial sound	Yes	No	No	No	No	No
Binaural Synchronization	Yes	Yes	No	No	No	No
Binaural Coordination	Via Streamer	Via Streamer	Via Streamer	No	No	No
Streamer/ConnectLine	Optional	Optional	Optional	No	No	No
Automatic Directionality	MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Single Band Adaptive Directionality; 3 modes: surround, split, full	MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	MultiBand Adaptive Directionality (4 Bands); 3 modes: surround, split, full	Single Band Adaptive Directionality; 3 modes: surround, split, full
Noise Management	TriState with Voicefinder	TriState with Voicefinder	Speech-weighted, modulation based	TriState with Voicefinder	TriState with Voicefinder	Speech-weighted, modulation based
Bandwidth	10kHz	8kHz	8kHz	10kHz	8kHz	8kHz
Adjustable Channels	10	8	6	10	8	6
DFC2 Feedback Management	Yes (binaural)	Yes (binaural)	Yes (binaural)	Yes	Yes	Yes
Voice Aligned Compression Fitting Rationales (Clarity2)	Yes	Yes	Yes	Yes	Yes	Yes
Identities	5	5	3	5	5	3
Front Focus	Yes	Yes	Yes	Yes	Yes	Yes
My Voice	Yes	No	No	No	No	No
Life Learning	Yes	No	No	Yes	No	No
Activity Analyzer	Yes	Yes	Yes	Yes	Yes	Yes
Number of Programs	Up to 4 via Streamer	Up to 4 via Streamer	Up to 4 via Streamer	1	1	1
AutoPhone and Auto T-Coil	Yes	Yes	Yes	Yes	Yes	Yes
Nanocoating	Yes	Yes	Yes	Yes	Yes	Yes
Battery size/hours	312 / (140 hours)	312 / (140 hours)	312 / (140 hours)	10a / (105 hours)	10a / (105 hours)	10a / (105 hours)
nEARcom	Yes	Yes	Yes	No	No	No